





ANNEK



VİZYON



VİZYON



ARTEK LİNEER



LORECA

AUDIOLOGY CABIN

REV	REVISIONS	DESCRIPTION	DATE

- ODYO KABIN (SESSİZ ODA)**
- Oda 120x120x220 cm dış ölçülerinde,
 - Panel Kalınlığı 15 cm,
 - Oda yüksekliği 220 cm,
 - 70x180cm çarpma kapı,
 - 15 cm plywood tabanlı,
 - Zamından 92cm yüksekliğe 70x50cm 6+12+6+12+6+12+6mm kalınlıkta 4 katlı cam pencere,
 - Pencere önüne 80x40cm katlanır masa konulacaktır.

THE INFORMATION AND THE ILLUSTRATIONS BELONG TO THE RIGHTS OF THE MANUFACTURER OR DISTRIBUTOR. THE INFORMATION IS PROVIDED AS IS AND WITHOUT WARRANTY OF ANY KIND. THE MANUFACTURER OR DISTRIBUTOR IS NOT RESPONSIBLE FOR ANY DAMAGE OR LOSS OF PROFITS OR REVENUE ARISING FROM THE USE OF THE INFORMATION. THE INFORMATION IS PROVIDED AS IS AND WITHOUT WARRANTY OF ANY KIND. THE MANUFACTURER OR DISTRIBUTOR IS NOT RESPONSIBLE FOR ANY DAMAGE OR LOSS OF PROFITS OR REVENUE ARISING FROM THE USE OF THE INFORMATION.

DOĞUŞ MEDİKAL TİC.VE SANLTD.ŞTİ.			
Description	Date	Dwg	
Audiology Cabin	18.08.2012	1/A	
Drawn by	Sevda	Şişek	
TG	1		1 OF 1